


PRODUCT **CATHETER PLUG**

- Raw Material:** Polypropylene (PP), DEHP Free, blue color
- Lenght:** 48 mm
- Diameter:** Universal, it fits any connector / catheter
- Tip:** Closed, progressively larger external diameter for perfect connection to catheters of any size
- Sterilization:** made of ETO + CO₂
- Registered Trademark:** Securdrain Vià!

DRAWING



CODE AND PACKAGE

CODE	DESCRIPTION	PACKAGE
11201102001*	Catheter Plug	100 pcs / box 3000 pcs / carton